

Routine production in a GMP environment of multiple ¹⁸F-radiopharmaceuticals on a disposable-based fully automated platform



Gameiro C.¹, Lambert B.², Sauvage C.³, Kozirowski J.⁴, Ilan Z.⁵, Schmitz A.⁶, Freifelder R.⁷, Ackermann U.⁷ and Tochon-Danguy H.J.⁷
¹IBA Chemin du Cyclotron, 3 Louvain-La-Neuve/BE ²IBA, 100 Executive Dr., Sterling, VA/ US ; ³IBA Pharma, Av.de l'Espérance, 1 Fleurus/BE; ⁴Herlev University Hospital Dep. Clinical Physiology, 54 P1 Herlev/DK; ⁵Aposense Ltd, 5-7 Odem St., Petach-Tikva/IL; ⁶University of Pennsylvania, Dep. of Radiology, Philadelphia PA/USA; ⁷Austin Health, Centre for Nuclear Medicine and PET, Heidelberg/AU

Introduction

PET modality is one of the most rapidly growing areas of medical imaging thanks to the availability of innumerable clinical centers with their own biomedical cyclotrons. To be able to cover the growing clinical demands a flexible, reliable platform (IBA Synthera) was developed for a (c)GMP environment. Well-known conventional chemistry steps were fully automated allowing the synthesis of multiple ¹⁸F-radiopharmaceuticals beyond FDG.



Synthera® Platform designed for a GMP environment

Materials/Methods


The platform consists of a synthesizer and a HPLC. The synthesizer employs a disposable system (IFP-integrated fluidic processor) where the entire synthesis takes place. The IFPs are named after conventional synthesis steps they are designed to perform: nucleophilic, chromatography, alkylation, distillation, and reformulation. One IFP only is needed for a single synthesis, but they can be connected in series for multi-step processes.

IFP™ cassette (Integrated Fluidic Processor)



Results

In this work, several ¹⁸F-labeled tracers were synthesized by nucleophilic fluorination (S_N2 and S_NAr) and synthon chemistry. The former can be performed in one step (¹⁸F-FP-DTBZ¹, ¹⁸F-fallypride²) or in two steps; fluorination followed by removal of protective groups (¹⁸F-FDG, ¹⁸F-FMISO³, ¹⁸F-FAZA², ¹⁸F-FHBG¹, ¹⁸F-Facetate⁴, ¹⁸F-F-A85380¹, ¹⁸F-ML-10, florbetaben). Synthons like ¹⁸F-FBM (fluorobromomethane), ¹⁸F-SFB⁴ (succinimidylfluorobenzoate) were also synthesized and they may be coupled to adequate precursor via acylation, alkylation or amide formation. For instance, ¹⁸F-FCH⁵ synthesis; ¹⁸F-FBM synthon produced in the first step N-alkylates the precursor resulting in ¹⁸F-FCH. This example shows the system is able to carry out multi-step process.



Synthera® & Synthera® HPLC Main page: fluidic pathway scheme, time counters ; signal recording window for real-time monitoring
 Report Page: report information can be customized, (from left to right)



Discussion/Conclusion

In most cases, the crude synthesis product required HPLC purification (*) while for the others cartridge purification (**) was sufficient. In every synthesis parameters were optimized with respect to precursor amount, reaction time, temperature and concentration. In the majority of the cases synthesis time was < 60 min. even when HPLC purification was included. Good synthesis yields and > 95 % radiochemical and chemical purity were obtained and were superior when compared to manual synthesis (yields at least doubled). Less radiation exposure, shorter synthesis time and stable yields are the other advantages versus manual synthesis. By simply adapting the recipe and by using adequate IFP the automated platform was able to consistently produce several tracers in high yields and suitable for human injection.